Fegura® Bond

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Instructions for use

English

Fegura® Bond is a bonding system with dual curing properties for making titanium and NP metals (on Co-Cr-Mo basis) adhere to prosthesis resins when realising removable dentures in the model casting technique.

The **Fegura® Bond** System is composed of:

- Primer
- Opaker powder – pink, dark pink, white, A3
- Liquid

Indication:

The **Fegura® Bond** Primer is used for conditioning surfaces of dental alloys in order to increase bonding and adherence of model casting alloys with denture resins on acrylic basis.

The **Fegura® Bond** Opaker powder and Liquid make it possible to adapt the colour of the metallic framework to the natural colour of gingiva and acrylic teeth.

Contra indications:

In case of a proven sensitizing or allergy against methyl methacrylate derivatives this bonding system should not be used, i.e. any processing of products should be replaced by other materials.

Instruction for use:

Preparing the casting model framework:

Cast the metallic framework as usual. For castings with a high fitting accuracy we recommend the **Feguramed model casting system**.

1. Blasting the framework:

Blast the framework with pure, clean aluminium oxide blasting materials (e.g. Alumix) with a grain size of 50 to 120 µm. Just blast the metallic frame with oil-free and steam-free compressed air for a short time.

2. Steam cleaning of the framework:

The blasted and dry metal framework should be steam-cleaned thoroughly with a steam jet or with pure acetone and dried carefully in the air. The surfaces to be coated must be perfectly dry and oil-free, otherwise, it could come to a premature wear of the metal/plastic bonding (fissure formation).

3. Application of the Fegura® Bond Primer:

Apply with the brush a layer of 100 to 300 µm (a well-covering brush coating) of **Fegura® Bond** Primer on the whole metal

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surface to be covered with dental resin after having roughened and cleaned it before (after use, immediately clean the brush with acetone!). Curing in the light curing unit should take place at least during 360 sec. with a Xenon flash bulb. We recommend the light curing units UniXS, HeraFlash (Heraeus), PLC Spectra (Schütz) and Spectra Link (Ivoclar) or a unit with the same radiation performance.

4. Mixing and application of the covering varnish:

Mix the opaker powder and the liquid in a clean small bowl in order to obtain a creamy consistency. First, take the required quantity of powder and wet it little by little with liquid and mix it until you obtain a consistency adapted for a brush application. Apply a thin, regular and well-covering layer with a clean brush (after use, immediately clean the brush with acetone!). The coloured coating layer must be cured at least during 360 sec. in a light curing unit, e.g. UniXS, HeraFlash (Heraeus), PLC Spectra (Schütz) or Spectra Link (Ivoclar). It is a **dual-curing** system, this means that polymerisation takes place at the same time under the effect of light and heat.

5. Dental resins for dentures:

In connection with the **Fegura® Bond** System we recommend following dental resins for dentures: Palapress (Heraeus), Futurapress (Schütz) and Degupress (Degudent).

Fegura® Bond Opaker powder and Liquid are foreseen for a single use only! A repeated use of the mixed polymer rests leads to a premature rupture of the metal/resin bonding!

Inhalation of acrylic resin vapours and dusts when processing or grinding is harmful to health! A direct contact of the skin, mucosa or eyes with not yet cured components can lead to a sensitising. Therefore, only process these products at working places provided with an exhaust system and avoid any direct contact with the skin, eyes and mucosa.

After use, immediately close the jar and store at 4 – 25 °C! Protect against light and heat! After the date of expiration – see the label on the packing – do not use any longer!

Package

REF

2180	Set	
2181	Primer	7 g
2182	Opaker powder pink	10 g
2182-2	Opaker powder dark pink	10 g
2183	Opaker powder white	5 g
2183-2	Opaker powder A3	5 g
2184	Liquid	10 ml